

# healthy skin rejuvenator step-by-step instructions

---



- 1.** Place the Puff on the Velcro Attachment Head — blue side up.



- 2.** Wet the Rejuvenating Puff, but don't drench it.



- 3.** Then snap the Attachment Head on to the System Applicator.



- 4.** Turn on the System Applicator and select your speed.



- 5.** Gently massage your face for 1-2 minutes per facial area, avoiding the eyes. Then rinse thoroughly.



- 6.** Remove and dispose of the Rejuvenating Puff and keep the Attachment Head. For best results, use Healthy Skin Rejuvenator 3 times a week.